










# High Quality, Durable Cryogenic & Laboratory Labels

perfect for proper labelling and accurate identification

- **Organise research processes** for better storage and identification
- **Assure marking legibility and labeling permanency** under harsh chemicals, extreme temperatures and extended storage
- **Brighten up the lab** with color-coded labels to differentiate processes and prevent errors

**40% OFF**

All Diversified Biotech Products

CRYOGENIC & LABORATORY LABELS	<b>CRYO-BABIES® &amp; CRYOTAGS®</b> 	<b>TOUGH-SPOTS®</b> 	<b>TOUGH-TAGS®</b> 
Ideal for	Cryogenic Vials Lg Vials/Tubes Microplates Racks & Boxes Bottles Microscope Slides	Cryogenic Vials Lg Vials/Tubes Microplates Racks & Boxes Bottles Microscope Slides	Tubes PCR Tubes Microscope Slides Glassware
Temperature Resistant			
Colours Available			
Adhere to	Plastic Glass Metal	Plastic Glass Metal	Plastic Glass Metal
<b>Marking Method</b>			
Handwritable	✓	✓	✓
Laser Printable	✓	✓	✓
Direct Thermal	✓		
Thermal Transfer	✓		✓

# FIND YOUR FIT

Sizing guide to match the perfect size and shape for your labelling needs.

## Sizing Guide (Actual size shown)

### TOUGH-TAGS®

24mm x 13mm  
For 0.5ml tubes

33mm x 13mm  
For 1.5-2.0ml tubes

Colors Available:

38mm x 19mm  
General Purpose Label

38mm x 6mm For Microplates

21mm x 7mm

Color Available:

## Sizing Guide (Actual size shown)

### TOUGH-SPOTS®

9.5mm Diameter  
For 0.5ml tubes

13mm Diameter  
For 1.5-2.0ml tubes

Colors Available:

4.8mm Diameter  
For 0.2ml tubes

11mm Diameter  
For 1.5-2.0ml tubes

19mm Diameter  
For 15-50ml tubes

25mm Diameter  
For 50ml tubes

Color Available:

## Sizing Guide (Actual size shown)

### CRYO-BABIES® AND CRYOTAGS®

38mm x 19mm  
General Purpose

33mm x 13mm  
For 1.5-2.0ml tubes

24mm x 13mm  
For 0.5ml tubes

43mm x 19mm  
CryoVials

38mm x 6mm (Microplates)

67mm x 25mm  
For boxes, tube racks & other containers

Colors Available:

76mm x 51mm  
For bottles, tube racks, boxes and  
other laboratory containers

64mm x 13mm  
For 0.5-2.0ml tubes

25mm x 25mm  
For Microscope  
Slides

44mm x 28mm  
For racks and boxes

Color Available:

Place your order with us today!

T. +65 6775 7318 E. [custcare@axilscientific.com](mailto:custcare@axilscientific.com) W. [www.axilscientific.com](http://www.axilscientific.com)