



- » Animal origin-free » Endotoxin tested » pH 7.1-7.5
- » Sterility tested for bacteria, fungus and mycoplasma



Medium Name	Product No.	Glucose	L-Glutamine	Phenol Red	HEPES	Sodium Pyruvate	Sodium Hydrogen Carbonate	Ribonucleoside	Deoxyribonucleoside	Pyridoxine Hydrochloride	Pyridoxal Hydrochloride	Non - Essential Amino Acids
DMEM	08456	1.0 g/l	○	○		○	○			○		
	08490	1.0 g/l				○	○			○		
	08457	4.5 g/l	○	○	25 mM		○			○		
	08458	4.5 g/l	○	○		○	○			○		
	08459	4.5 g/l	○	○			○			○		
	08489	4.5 g/l					○			○		
	08488	4.5 g/l		○			○			○		
	11584	4.5 g/l		○			○			○		
	11585	4.5 g/l		○		25 mM		○		○		
DMEM/Ham's F-12	09891	No	○	○			○			○		
	08460	3.15 g/l	○	○	15 mM	○	○			○		○
	05177	3.15 g/l	○		15 mM	○	○			○		○
	11581	3.15 g/l	○	○		○	○			○		○
	11582	3.15 g/l	○			○	○			○		○
	11583	3.15 g/l		○	15 mM	○	○			○		○
Ham's F-12	09893	No	○	○		○	○			○		○
	17458	1.8 g/l	○	○		○	○			○		○
IMDM	11506	4.5 g/l	○	○	25 mM	○	○				○	○
MEM	21442	1.0 g/l	○	○			○				○	
	21443	1.0 g/l	○	○			○				○	○
	09848	No	○	○			○				○	○
α-MEM	21444	1.0 g/l	○	○		○	○	○	○		○	○
	21445	1.0 g/l	○	○		○	○				○	○
RPMI1640	30263	2.0 g/l	○	○	25 mM							
	30264	2.0 g/l	○	○			○			○		
	06261	2.0 g/l	○				○			○		
	05176	2.0 g/l		○			○			○		
	09892	No		○	○			○		○		

DMEM/Ham's F-12

"So far after trying out the 08460-95 Media on the cells, the cells are well maintained. The ordering and delivery of the product are efficient and flexible as well. Good job Axil Scientific team!"

- Research Officer, Biopolis

RPMI 1640

"Cited in a lot of literature. Superior cell growth and efficient delivery system"

- Research Assistant, NUS

DMEM / RPMI 1640

"Superior cell growth for breast cancer cells 231 mcT7"

- Researcher, NUS

DMEM

"Nacalai has been a very new brand to me and I wasn't really sure if I should change my current brand. However after using the sample bottle of DMEM provided by Axil Scientific, my cells continued to grow healthy."

- Research Assistant, Biopolis